


11113

Kit for metal detection, 55 mm, Blue



Use this test kit to check the detectability of Vikan metal-detectable plastic fragments as part of your HACCP plan.

Technical Data

Item Number	11113
Material	Polypropylene Metal & X-Ray detectable additive
Complies with (EC) 1935/2004 on food contact materials ¹	Yes
Complies with EU Regulation 2023/2006/EC of Good Manufacturing Practice	Yes
Complies with FDA Regulationl CFR 21 ¹	Yes
Complies with UK 2019 No. 704 on food contact materials	Yes
Complies with REACH Regulation (EC) No. 1907/2006	Yes
Complies with California Proposition 65	Yes
Complies with Halal and Kosher	Yes
PFAS, Phthalates and BPA intentionally added	No
Box Quantity	5 Pcs.
Quantity per Pallet (80 x 120 x 180-200 cm)	50000 Pcs.
Quantity Per Layer (Pallet)	0 Pcs.
Box Length	0 mm
Box Width	0 mm
Box Height	0 mm
Length/Depth	55 mm
Width	45 mm
Height	8 mm
Net Weight	0.004 kg
Weight bag (Recycling Symbol "4") LDPE ≥ PCR 30%	0.005 kg
Weight cardboard (Recycling symbol "20" PAP)	0.0058 kg
Tare total	0.0108 kg
Gross Weight	0.01 kg
Cubic metre	2E-05 M3
Recommended sterilisation temperature (Autoclave)	121 °C
Max. usage temperature (food contact)	100 °C
Max usage temperature (non food contact)	100 °C
Min. usage temperature ³	-20 °C
Min. pH-value in usage concentration	2 pH
Max. pH-value in Usage Concentration	10.5 pH
Metal Detectable	Yes
GTIN-13 Number	5705022033498
GTIN-14 Number (Box quantity)	15705028033505
Customs Tariff No.	39249000

New equipment should be cleaned, disinfected, sterilised and any labels removed, as appropriate to its intended use, e.g. high risk vs. low risk food production areas, general hospital areas vs. intensive care units, before use.

1. See Declaration of Compliance for further details on food contact
3. Do not store the product below 0° Celsius.